採尿方法 How to take urine

別紙4 Separate Sheet No.4

1) 採尿カップに尿を取る
Take urine in a
test cup.



2) 採尿容器の腹を 押さえ、採尿カップの 尿を吸い上げる



Pushing the side of the bottle, pump up urine from the cup by a bottle.



3) 線のところまで 吸い上げ、容器の ふたを閉める

Fill the bottle to a marker line, and screw the cap back on.



4) 提出用袋に 所属と氏名を記入 し提出

Write your affiliation and name on the bag.

